

## TECHNICAL DATA SHEET

# SteriCleen Neutral Detergent 5L

**Neutral Medical Instrument & Equipment Cleaner**  
**Australian Medical Device Class 1**  
**ARTG No. 332455**

**CODE: 9220/42**

### DESCRIPTION:

SteriCleen Neutral Detergent is a highly active neutral medical device detergent specifically designed for the removal of hospital and surgical soils via both machine and manual cleaning operations. SteriCleen Neutral Detergent is a concentrated formula which is low-foaming, non-staining and effective in both hard and soft water.

### BENEFITS:

- Concentrated blend with neutral pH– safe for use on most instruments and devices
- Biodegradable surfactant system with recyclable packaging
- Contains no ammonia or phosphates
- Non-toxic, non-abrasive, low foaming, free rinsing
- Contains corrosion inhibitors for the prevention of corrosion of medical instruments during the cleaning process
- Can be used in Ultrasonic & Manual cleaning of instruments prior to disinfecting or sterilisation procedures
- Complies with AS/NZS4187 and AS/NZS4815

### INTENDED USAGE AND DIRECTIONS TO USE:

- Use as a dedicated cleaner in conjunction with Ultrasonic cleaning machines
- Manual cleaning of Medical and Dental instruments prior to the disinfection or sterilisation process
- Manual cleaning of theatre floors, tables, benches as well as laboratory glassware
- Manual cleaning of large and small equipment

### DIRECTIONS FOR USE:

1. Manufacturer's instructions should be strictly followed as advised by the instrument provider.
2. Prepare a fresh ready to use solution daily.
3. Dilution prior to use:
  - Heavy soiling: Dilute 10 mL per Litre of water
  - Medium to light soiling: Dilute 5 mL per Litre of water
 (Discard diluted solutions within 24 hours of preparation)

Note: Instrument cleaners are not disinfectants, and all instruments should be cleaned, rinsed and drained before immersion in an instrument grade DISINFECTANT, or prior to autoclaving.



## TECHNICAL DATA SHEET

### CHEMICAL COMPOSITION AND PROPERTIES:

Colour concentrate:	Yellow
Colour dilute solution:	Colourless (5 mL per litre)
Odour (Concentrate):	No odour
Surfactant system:	Non-ionic
pH concentrate:	7.0-9.0
pH dilution (5 mL per litre):	7.5 ± 0.5
Solubility:	Completely soluble in water
Rinsibility:	Residue free
Specific gravity:	1.00 ± 0.03g/ml
Corrosion inhibitor:	Proprietary blend
Shelf life:	3 years

**STORAGE CONDITIONS:** Store below 30°C.

**SAFETY DATA:** Information is listed in the Material Safety Data Sheet.

**PACK SIZE:** 5 Litres

**WARRANTY** – All statements, information and data presented herein are believed to be accurate and reliable but are not to be taken as a guarantee, expressed or implied, for which seller assumes legal responsibility and they are offered solely for your consideration, investigation and verification. Statements or suggestions concerning possible use of this product are made without representation or warranty that any such use is free of patent infringement and are not recommendations to infringe on any patent. Created 30 November 2021 Date Printed 15/12/2021 1:52 PM